


cuesee

Spotlight

Oct 2024



Eurotrol

eurotrol.com



Quality Controls for In Vitro Diagnostics

For 30+ years Eurotrol has been entrusted by customers worldwide to develop, manufacture and distribute CueSee® solutions for In Vitro Diagnostics (IVD).

Typical customers are


- Hospital laboratories
- Point-of-care sites
- Proficiency testing organizations
- Standardization bodies
- Diagnostic companies

Products have been registered in more than 100 countries.

We do offer CueSee® customized to meet your specific needs.



Blood Gases



Cardiac Markers



Diabetes Care



Hemoglobin



Coagulation

On the following pages you will find spotlights of some latest developments in CueSee®.

Eurotrol


For a more detailed overview of our products and developments, please refer to: eurotrol.com/products





Your Companion
to Validate
CueSee® Results

There is no more reliable reference for the expected results of CueSee® than peer data comparison. Eurotrol offers a comprehensive and user friendly solution to assist you with this reference for all your CueSee® testing needs. There are endless possibilities for statistical and graphical presentations based on your own acceptance criteria. You will soon discover that it is the perfect companion for your CueSee®.



Sign up for CueSee® Compare NOW!

cuseecompare.com

For any assistance contact our team
supportus@eurotrol.com | 502-501-1180 Ext. 1

Benefits


- Comparison of your results with peers
- Simple result entry manually or import
- Multiple reports with **your acceptance criteria**
- Instant access to all historical reports
- Customizable to your testing program
- Satisfy regulatory requirements
- Continuous training and support
- **FREE** with all Eurotrol products



CueSee® Ringe is an innovative packaging concept that will modernize CueSee® for blood gases on all iSTAT BGEM cartridges. A special glass ampoule sealed inside the device is opened with a simple turning motion eliminating the risk of cut fingers and glass scattering. After activation it functions similar to a standard syringe. The liquid inside the ampoule is not contaminated by ambient air, which reduces pre-analytical errors and provides a long open vial stability. CueSee® Ringe will soon be the new standard in CueSee® of blood gases.



CueSee VeriSTAT® QC Level | 1 | 2 | 3 |



CueSee VeriSTAT® LV Level | 1 | 2 | 3 | 4 | 5 |

Benefits

- Safe and easy sampling
- No shaking required before use
- Long shelf-life of 4 years at 2-8 °C
- 30 days storage at room temperature before use
- All analytes in one vial including Hematocrit (Hct)
- **Blue dye** to enhance sample visibility
- Products are offered as a QC and LV bundle
- The best precision for gases
- **Open vial stability of 15 minutes!**
- Multiple tests with just ONE CueSee® Ringe

Ideal for routine CueSee®, method comparisons, reportable range validation, linearity studies and competency testing in combination with **FREE** access to CueSee Compare.





Try it and CueSee®
will never be the
same again

The only CueSee® for CO-oximeters covering the whole reportable range with real hemoglobin derivatives like patient samples. Packaged in the ACU-Drop II, a two chambered device, which separates the FO₂Hb and FCOHb from FMetHb during storage, provides a stable formulation. This separation of fractions also means multiple concentrations of all analytes including FMetHb are available. Once ready to use, simply push a button to mix the fractions and use the built-in dropper bottle or attach a syringe to test.



CueSee CO-OX QC Level | 1 | 2 | 3 |



CueSee CO-OX VM Level | 1 | 2 | 3 | 4 | 5 | 6 |

Benefits

- Containing measurable ranges for tHb, FO₂Hb, FCOHb and **FMetHb (Methemoglobin)** for the whole reportable range
- Safe and easy sampling
- Open vial stability of 15 minutes
- Proven superior performance on CO-oximeters
- Products are offered as QC and VM bundle
- Rated no. 1 by Proficiency Testing (PT) organizations

Ideal for routine CueSee®, method comparisons, reportable range validation, linearity studies and competency testing in combination with **FREE** access to CueSee Compare.





We made the
impossible possible

Eurotrol has revolutionized CueSee® for coagulation by developing a unique plasma with red blood cells for reporting PT, INR, ACT and aPTT. Packaged in the innovative CueSee® Mix, a dual-chambered device, that keeps the lyophilized plasma with red blood cells separated from the reconstitution solution until activation. When ready to test a simple twisting motion allows the two components to mix. After mixing, the CueSee® material is easy and safe to test by use of the integrated dropper or by attaching a syringe.



CueSee® Coag QC Level | 1 | 2 | 3 |



CueSee® Coag VM Level | 1 | 2 | 3 | 4 | 5 |

Benefits

- Safe and easy sampling
- Covering the clinical relevant ranges for PT, INR, ACT and aPTT in **ONE** vial
- Plasma with red blood cells that mimic patient samples
- Products are offered as a QC and VM bundle
- Compatible with all common POCT devices

Ideal for routine CueSee®, method comparisons, reportable range validation, linearity studies and competency testing in combination with **FREE** access to CueSee® Compare.



Meet us

Medica 2024

11-14 Nov 2024 Düsseldorf Germany

With over 5,300 exhibitors from almost 70 countries and 83,000 visitors MEDICA in Düsseldorf is one of the largest medical B2B trade fairs in the world. A wide range of innovative products and services from the fields of medical imaging, laboratory technology, diagnostics, health IT, mobile health as well as physiotherapy / orthopaedic technology and medical consumables are presented here.

AARC Congress

20-23 Nov 2024 Orlando, FL USA

The AARC International Congress is the premier networking and educational event for respiratory therapists from across the country and around the world. Congress offers cutting-edge lectures by the top leaders in the industry, a robust exhibit hall showcase the latest in industry products and services, and, of course, ample networking.


COLA

6-8 May 2025 Chandler, AZ USA

The educational sessions at the Laboratory Enrichment Forum provide practical tools to use in your laboratory as you support other medical professionals and the communities you serve together. Industry leaders and experts will discuss laboratory leadership, quality improvement, emerging trends and new technologies, leaving you better equipped to thrive in an unpredictable landscape. Continuing education credits are available for physicians and laboratorians.



Follow us



in



f



x

cuesee

Find us

Headquarters

Eurotrol B.V.
Keplerlaan 20
6716 BS Ede
The Netherlands
office@eurotrol.com
T +31 318 695777

The Americas


Eurotrol, Inc.
850 N Black Branch Rd
Elizabethtown, KY 42701
United States
officeus@eurotrol.com
+1 502 501 1180



eurotrol.com

cuesee

Spotlight



Eurotrol